

# Panther has new prey.



**A new standard: the ideal combination of performance and automation for HIV-1 viral load monitoring.**

Excellent sensitivity and high precision.<sup>1</sup>

Performance	
<b>LoD</b>	12 copies/mL (3rd WHO International Standard)
<b>LLoQ</b>	30 copies/mL (Reaction volume: 0.5 mL plasma)
<b>Linear range</b>	30 to 10 million copies/mL
<b>Accuracy</b>	External quality assessment panels used to verify accuracy across known standards—AcroMetrix, Qnostics, SeraCare, CAP, QCMD

Mean concentration (Log copies/mL)	Total	
	SD	CV %
1.80	0.16	9.10
2.37	0.17	7.34
2.47*	0.15	6.15
2.95	0.15	5.20
3.80	0.10	2.53
4.93	0.08	1.56
5.69	0.09	1.58
6.71	0.08	1.13


\*Panel member was diluted 1:3 just prior to run. Reported result (shown here) has dilution factor of 3 applied automatically by the Panther® system.

Excellent precision with the Aptima HIV-1 Quant assay provides confidence that increases in viral load are due to clinical changes with the patient versus assay variation.<sup>2†</sup>

†Only the data with respect to the viral load quantitation is applicable to the Aptima HIV-1 Quant assay. This product is marketed as the Aptima HIV-1 Quant Dx assay outside of the US, because it is also intended for diagnosis of HIV-1 outside of the US.

**Demand more from your HIV-1 assay.**

Product design	
<b>Intended use</b>	▶ HIV-1 viral load monitoring
<b>Technology</b>	▶ Real-time transcription-mediated amplification (TMA)
<b>Target region</b>	▶ pol/LTR
<b>Groups</b>	▶ M, N, O
<b>Sample types</b>	▶ Quantitative: Plasma
<b>Sample input volume</b>	▶ Primary tube (PPT, ACD, EDTA): 1.2 mL ▶ Secondary tube: 700 µL ▶ 240 µL with automated dilution factor



**Aptima® HIV-1**  
Quant Assay

# The Aptima<sup>®</sup> HIV-1 Quant assay on the Panther<sup>®</sup> system

## The ideal combination of performance and superior automation for viral load monitoring of HIV-1.

Sample to result within a single integrated instrument.

Key automation characteristics	
<b>Random access</b>	<ul style="list-style-type: none"> <li>▶ No more batching; load samples with different test orders as they arrive</li> <li>▶ Run HIV with other assays at the same time</li> </ul>
<b>Primary tube processing</b>	<ul style="list-style-type: none"> <li>▶ No need for aliquoting or manual sample transfer</li> <li>▶ Tube flexibility: PPT, ACD and EDTA tubes validated</li> </ul>
<b>Flexible sample and reagent loading</b>	<ul style="list-style-type: none"> <li>▶ Run multiple tests from same sample</li> <li>▶ No manual sample prep or barcode clips</li> <li>▶ Automated barcode scanning of tubes allows positive sample identification</li> </ul>
<b>Rapid turnaround time with stat result option</b>	<ul style="list-style-type: none"> <li>▶ First results in 2 hours 41 minutes</li> <li>▶ Stat results option: ability to prioritize results</li> </ul>
<b>Automated QC analysis</b>	<ul style="list-style-type: none"> <li>▶ Levey-Jennings plots to track and trend controls</li> </ul>
<b>Low sample volume option (240 µL) with automated dilution factor</b>	<ul style="list-style-type: none"> <li>▶ Specified dilution factor automatically applied versus manual calculations</li> </ul>

## Ordering information.

Item	Quantity	Catalog #
Aptima HIV-1 Quant Assay Kit (1 assay box, 1 calibrator kit, 1 control kit)	100 tests	PRD-03565
Aptima HIV-1 Quant Controls Kit	5 of each control (low, negative, high)	PRD-03567
Aptima HIV-1 Quant Calibrator Kit	5 calibrators	PRD-03566

For more information contact your Hologic representative or visit [USAptimaVirology.com](http://USAptimaVirology.com).

### References:

1. Aptima HIV-1 Quant Assay [package insert] #AW-13242-001, San Diego, CA; Hologic, Inc., 2016. 2. Nair SV, Kim HC, Fortunko J, et al. Aptima HIV-1 Quant Dx—a fully automated assay for both diagnosis and quantification of HIV-1. *J Clin Virol*. 2016;77:46-54.

SS-00425-001 Rev. 001 ©2016 Hologic, Inc. All rights reserved. Hologic, The Science of Sure, Aptima, Panther, and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to [diagnostic.solutions@hologic.com](mailto:diagnostic.solutions@hologic.com).

**Aptima<sup>®</sup> HIV-1**  
 Quant Assay